




# X-RAY DETECTABLE THREAD



It can be **weaved or sewed** (or welded by using ultrasonic waves or heater) into gauze, non-woven or operation towel.

Whether the gauze or the operation towel is left in the body of a patient can be found through X-ray examination after operations, and medical negligence caused by the gauze or the operation towel left in the body of the patient is **avoided**.

Copyright by LOYAL. This brochure contains photos and technical data for information only, without contractual engagement. LOYAL reserves the right to introduce any modifications in the design and/or equipment.



## WUHAN LOYAL MACHINERY CO., LTD

© 5B05, Decheng Center (DECH CENTER), Zhongbei Road No. 124,  
Wuchang District, 430061, Wuhan City, Hubei Province, China  
Mobile: 0086 1387 1887 978 (WhatsApp / Wechat) Tel / Fax: 0086 27 8566 3380  
E-mail: loyal@med-machinery.cn michelle\_6237@hotmail.com  
www.med-machinery.cn



# X-RAY DETECTABLE THREAD



**Loyal** accept customized different diameter range from **0.55-0.68mm**.

Accept customized different melt point from **135-175°C**, in case to meet different medical dressing.

Our relevant machines can do customize for your on-line perm X-ray detectable thread on dressing request, please feel free to contact us for more information.

Copyright by LOYAL. This brochure contains photos and technical data for information only, without contractual engagement. LOYAL reserves the right to introduce any modifications in the design and/or equipment.



**WUHAN LOYAL MACHINERY CO., LTD**

© 5B05, Decheng Center (DECH CENTER), Zhongbei Road No. 124,  
Wuchang District, 430061, Wuhan City, Hubei Province, China  
Mobile: 0086 1387 1887 978 (WhatsApp / Wechat) Tel / Fax: 0086 27 8566 3380  
E-mail: loyal@med-machinery.cn michelle\_6237@hotmail.com  
www.med-machinery.cn