




# SERIE SS-FQ

*Gauze Slitting and Re-winding Machine*



Fully automatic slitting and re-winding machine, working from original gauze rolls (from loom or after bleaching process).



Specially designed to obtain pre-cutted gauze rolls to be used in the feeding of our automatic gauze swab folding machine (SERIE SS-ZDG) for the production of swabs ( compresses ).

Copyright by LOYAL. This brochure contains photos and technical data for information only, without contractual engagement. LOYAL reserves the right to introduce any modifications in the design and/or equipment.



**WUHAN LOYAL MACHINERY CO., LTD**

© 5B05, Decheng Center (DECH CENTER), Zhongbei Road No. 124,  
 Wuchang District, 430061, Wuhan City, Hubei Province, China  
 Mobile: 0086 1387 1887 978 (WhatsApp / Wechat) Tel / Fax: 0086 27 8566 3380  
 E-mail: [loyal@med-machinery.cn](mailto:loyal@med-machinery.cn) [michelle\\_6237@hotmail.com](mailto:michelle_6237@hotmail.com)  
[www.med-machinery.cn](http://www.med-machinery.cn)



# SERIE SS-FQ

## *Gauze Slitting and Re-winding Machine*


### TECHNICAL SPECIFICATIONS



The raw material holder has regulating function which can adjust the position of gauze jumbo roll position within range of  $\pm 2.5\text{cm}$ .

After gauze jumbo roll loading on machine, user can easily adjust the position without to lift gauze jumbo roll (very heavy).

Hauling and pushing reasonable combined together, maximum reduce the deformation when gauze been drive and movement over.



During winding, the tension control device works to make it suitable to cut raw gauze with too many joints and unbalanced tension.

Copyright by LOYAL. This brochure contains photos and technical data for information only, without contractual engagement. LOYAL reserves the right to introduce any modifications in the design and/or equipment.





**WUHAN LOYAL MACHINERY CO., LTD**

© 5B05, Decheng Center (DECH CENTER), Zhongbei Road No. 124,  
Wuchang District, 430061, Wuhan City, Hubei Province, China  
Mobile: 0086 1387 1887 978 (WhatsApp / Wechat) Tel / Fax: 0086 27 8566 3380  
E-mail: loyal@med-machinery.cn michelle\_6237@hotmail.com  
www.med-machinery.cn



# SERIE SS-FQ

## *Gauze Slitting and Re-winding Machine*




Both Pneumatic cutter or Electrical cutter can be customized on machine, user can according production request to choose the suitable cutter for machine.



By adjusting the position of the blade, user can easily control the cutting width.

According to production requirements, the online melting X-ray detectable yarn function can be added, it makes the production process more diversified(customization option).



Copyright by LOYAL. This brochure contains photos and technical data for information only, without contractual engagement. LOYAL reserves the right to introduce any modifications in the design and/or equipment.



**WUHAN LOYAL MACHINERY CO., LTD**


© 5B05, Decheng Center (DECH CENTER), Zhongbei Road No. 124,  
Wuchang District, 430061, Wuhan City, Hubei Province, China  
Mobile: 0086 1387 1887 978 (WhatsApp / Wechat) Tel / Fax: 0086 27 8566 3380  
E-mail: loyal@med-machinery.cn michelle\_6237@hotmail.com  
www.med-machinery.cn





# SERIE SS-FQ

## *Gauze Slitting and Re-winding Machine*



Frequency converter controls the speed, make constant tension and linear velocity, precise speed regulating, stable running, and accurate counting.

### Technical Parameter

Jumbo roll max. diameter	87cm
Max. working width	90-240cm (Accept customize)
Cutting speed	0-180m/min
Hauling motor power	4 KW
Pushing motor power	1.1 KW
Frequency converter	3.7KW; 5.5KW

Copyright by LOYAL. This brochure contains photos and technical data for information only, without contractual engagement. LOYAL reserves the right to introduce any modifications in the design and/or equipment.



**WUHAN LOYAL MACHINERY CO., LTD**

© 5B05, Decheng Center (DECH CENTER), Zhongbei Road No. 124,  
 Wuchang District, 430061, Wuhan City, Hubei Province, China  
 Mobile: 0086 1387 1887 978 (WhatsApp / Wechat) Tel / Fax: 0086 27 8566 3380  
 E-mail: loyal@med-machinery.cn michelle\_6237@hotmail.com  
 www.med-machinery.cn