

SERIE SS-DL

Single Wheel Folding and Rewinding Machine



Series single wheel folding and rewinding machine folds the gauze lengthwise into several plies and rewinds for the preparation of swabs.

It is easily to switch machine folding mode to achieve different pre-folding plies request.

Tension controlled by magnetic powder clutch, folding tension stable.

Frequency timing system makes machine operate easy and efficient.



 SB05, Decheng Center (DECH CENTER), Zhongbei Road No. 124, Wuchang District, 430061, Wuhan City, Hubei Province, China
Mobile: 0086 1387 1887 978 (WhatsApp / Wechat) Tel / Fax: 0086 27 8566 3380 E-mail: loyal@med-machinery.cn michelle_6237@hotmail.com www.med-machinery.cn



SERIE SS-DL

Single Wheel Folding and Rewinding Machine

TECHNICAL SPECIFICATIONS



made by this machine.

This machine is indispensable for use SS-ZD gauze swab folding machine (unfolded edge).




 5B05, Decheng Center (DECH CENTER), Zhongbei Road No. 124, Wuchang District, 430061, Wuhan City, Hubei Province, China
Mobile: 0086 1387 1887 978 (WhatsApp / Wechat) Tel / Fax: 0086 27 8566 3380 E-mail: loyal@med-machinery.cn michelle_6237@hotmail.com www.med-machinery.cn




Copyright by LOYAL. This broch

SERIE SS-DL

Single Wheel Folding and Rewinding Machine





On-line perm X-RAY detectable thread function can be customized.

Technical Parameter

Folding speed	20-150 meters/min
Working width	15cm, 22.5cm, 30cm (accepted customized)
Folding plies	2ply and 3 ply
Max. rewinding diameter	85cm

WUHAN LOYAL MACHINERY CO., LTD

contains photos and technical data for information only, without contractual engagement. LOYAL reserves the right to introduce any modifications in the design and/or equipment

 SB05, Decheng Center (DECH CENTER), Zhongbei Road No. 124, Wuchang District, 430061, Wuhan City, Hubei Province, China
Mobile: 0086 1387 1887 978 (WhatsApp / Wechat) Tel / Fax: 0086 27 8566 3380 E-mail: loyal@med-machinery.cn michelle_6237@hotmail.com www.med-machinery.cn